Upute korisnicima laboratorijskih usluga vezano za pretrage koje se izvode u dijagnostičkim laboratorijima iz područja citologije, imunofenotipizacije, citogenetike, molekularne dijagnostike i tipizacije tkiva

Toni Božac¹, Jozefina Palić¹, Iva Semren¹, Monika Kolundžić¹, Ivona Arić Zrna², Antonija Babić¹, Ivana Franić Šimić¹, Branimir Gizdić³, Koraljka Gjadrov Kuveždić⁴, Zorana Grubić⁵, Mirjana Marijana Kardum Paro⁶, Josip Knežević¹, Lejla Kurić⁶, Saška Marcić⁷, Sanja Mrsić Davidović¹, Marijana Nehutni¹, Sunčica Ries⁴, Mirjana Suver Stević⁷, Zoran Šiftar⁶, Katarina Štingl Janković⁵, Željka Švabek Tkalić⁶, Klara Zubčić², Marko Žarak², Renata Žunec⁵, Margareta Radić Antolčić¹

¹ Klinički zavod za laboratorijsku dijagnostiku, KBC Zagreb, Zagreb, Hrvatska

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ZA SVA DODATNA PITANJA KONTAKTIRAJTE NADLEŽNI SPECIJALIZIRANI LABORATORIJ!

Literatura

- Lippi G, Betsou F, Cadamuro J, et al. Preanalytical challenges - time for solutions. *Clin Chem Lab Med*. 2019;57(7):974-981. doi:10.1515/cclm-2018-1334
- Apperley JF, Milojkovic D, Cross NCP, et al. 2025 European LeukemiaNet recommendations for the management of chronic myeloid leukemia. *Leukemia*. 2025;39(8):1797-1813. doi:10.1038/s41375-025-02664-w
- Cloos J, Valk PJM, Thiede C, et al. 2025 update on MRD in acute myeloid leukemia: a consensus document from the ELN-DAVID MRD Working Party. *Blood*. 2026;147(11):1147-1167. doi:10.1182/blood.2025031480
- Gökbuğut N, Boissel N, Chiaretti S, et al. Diagnosis, prognostic factors, and assessment of ALL in adults: 2024 ELN recommendations from a European expert panel. *Blood*. 2024;143(19):1891-1902. doi:10.1182/blood.2023020794
- Thiele J, Kvasnicka HM, Orazi A, et al. The international consensus classification of myeloid neoplasms and acute Leukemias: myeloproliferative neoplasms. *Am J Hematol*. 2023;98(1):166-179. doi:10.1002/ajh.26751
- Agathangelidis A, Chatzidimitriou A, Chatzikonstantinou T, et al. Immunoglobulin gene sequence analysis in chronic lymphocytic leukemia: the 2022 update of the recommendations by ERIC, the European Research Initiative on CLL. *Leukemia*. 2022;36(8):1961-1968. doi:10.1038/s41375-022-01604-2
- Malcikova, J., Pavlova, S., Baliakas, P. et al. ERIC recommendations for TP53 mutation analysis in chronic lymphocytic leukemia—2024 update. *Leukemia* 38, 1455–1468 (2024). <https://doi.org/10.1038/s41375-024-02267-x>
- ELN - WP10 - Consensual recommendations on preanalytical precautions for the immunophenotyping of leukemia and immunoproliferative disorders. 2005 Oct 31. www.leukemia-net.org
- Davis BH, Dasgupta A, Kussick S, Han JY, Estrellado A; ICSH/ICCS Working Group. Validation of cell-based fluorescence assays: practice guidelines from the ICSH and ICCS - part II - preanalytical issues. *Cytometry B Clin Cytom*. 2013 Sep-Oct;84(5):286-90. doi: 10.1002/cyto.b.21105.
- Flores-Montero J, Sanoja-Flores L, Paiva B, Puig N, García-Sánchez O, Böttcher S et al. Next Generation Flow for highly sensitive and standardized detection of minimal residual disease in multiple myeloma. *Leukemia* 2017; 31: 2094-2103. doi:10.1038/leu.2017.29
- Hallek M, Cheson BD, Catovsky D, et al. iwCLL guidelines for diagnosis, indications for treatment, response assessment, and supportive management of CLL. *Blood*. 2018 Jun 21;131(25):2745-2760. doi: 10.1182/blood-2017-09-806398.
- Zwaan C. M., Tasian S. K., Aplenc R., et al. Diagnosis and Management of AML in Pediatric Patients: Consensus Recommendations from an International Expert Panel. *Blood*. 2025;blood.2024027904. doi:10.1182/blood.2024027904
- Hastings R, Howell R, Betts D, et al. Guidelines and quality assurance for acquired cytogenetics: A common European framework for quality assessment for banded chromosome studies and molecular cytogenetic investigations of acquired abnormalities. E.C.A. Permanent Working Group for Cytogenetics and Society; 2013.
- Mikhail FM, Heerema NA, Rao KW, et al. Section E6.1-6.4 of the ACMG technical standards and guidelines: chromosome studies of neoplastic blood and bone marrow-acquired chromosomal abnormalities. *Genet Med*. 2016 Jun;18(6)
- Shomali W, Gotlib J. World Health Organization-defined eosinophilic disorders: 2022 update on diagnosis, risk stratification, and management. *Am J Hematol*. 2022 Jan 1;97(1):129-148. doi: 10.1002/ajh.26352.